

# HeatSensor

HEATSENSOR® DUO - OCTO are incubators ensuring all Unisensor diagnostic tests to receive a fixed incubation temperature no matter what the ambient / operating temperature.

Both incubators are a **One Step Test Procedure**.

## KEY-POINTS

- » Automatic release of dipstick.
- » Compatibility with all test kits from UNISENSOR.
- » DUO is for maximum 2 tests at a time.
- » OCTO is for maximum 8 tests at a time.
- » Self-timing with light and audio signal.
- » Countdown timer display.
- » Product's specific program included.
- » Easy to use.



### SPECIFIC PROGRAM

Select the right program and place the dipstick.



### AUTOMATIC

The dipstick falls down automatically into the microwell after the first incubation time.

## PRODUCTS REFERENCES

**APP032**  
HEATSENSOR DUO

**APP052**  
HEATSENSOR OCTO



Unisensor s.a.

Liège Science Park  
Rue Louis Plescia, 8  
B-4102 Ougrée (Seraing)  
BELGIUM

[www.unisensor.be](http://www.unisensor.be)

Phone +32 4 252 66 02  
[info@unisensor.be](mailto:info@unisensor.be)

Reseller