

ReadSensor

READSENSOR® 2 is a handheld instrument designed to read results from all UNISENSOR test strips.

It quantifies, interprets, and records test strip results.

Simple to operate, READSENSOR® 2 provides rapidity, objectivity, convenience and the ability to store test data on your personal computer, or to print them with an embedded printer.

KEY-POINTS

- Clear rapid and accurate interpretation of results.
- Great repeatability of results.
- Results can be displayed, printed, stored or transferred to USB upon request.
- Compatibility with all test kits from UNISENSOR.
- Wide touchscreen.
- Operator Name and Product Lot number entry for increased traceability.
- Auto-calibration with external Reference dipstick.
- Automatic, motorized strip drawer.



EMBEDDED PRINTER

Results can be printed with the use of the embedded printer. A small ticket will ensure a good traceability of all results.



FASTER RESULTS

Results display in 5 seconds: Insert dipstick, select program and start reading.

PRODUCT REFERENCE

APP088
READSENSOR 2



Unisensor S.A.
Liège Science Park
Rue Louis Plescia, 8
B-4102 Seraing (Ougrée)
BELGIUM

Phone +32 4 252 66 02
info@unisensor.be
www.unisensor.be

Reseller